



SUTURES



There's sharp, and then there's TBS Graphite Black™ sharp.

TBS GRAPHITE BLACK™ NEEDLE

PRECISION REVERSE CUTTING

Setting you up for ultra sharp, anti-glare, multi-pass performance by default.

7-STAGE SHARPENING AND FINISHING PROCESS

Proprietary black oxidization process for maximum depth-perception and visibility, ultra-refined to puncture tissue effortlessly.

SILICONE COATING

Reduces tissue drag for buttery smooth punctures, pass after pass.

Precision Reverse Cutting



Introducing TBS Sutures - where reliability meets performance.

Featuring the ultra-sharp, high-contrast, TBS Graphite Black™ Needle for clean, consistent closures.

Now shipping with select models, these needles pierce effortlessly, reduce tissue drag, and deliver unmatched visibility during suturing.

Packaged in our tangle-free racetrack packaging with a full portfolio of absorbable and non-absorbable configurations, TBS Sutures make every stitch count.

POLYGLYCOLIC ACID

ABSORBABLE (PGA)

BRAIDED AND COATED

\$219.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-PGA3024GB	19mm	3/8	Black	3-0	24"	Violet
S-PGA4024GB	19mm	3/8	Black	4-0	24"	Violet

POLYGLACTIN 910

ABSORBABLE (PGLA)

BRAIDED AND COATED

\$219.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-PGLA3024GB	19mm	3/8	Black	3-0	24"	Violet
S-PGLA4024GB	19mm	3/8	Black	4-0	24"	Violet

PTFE (PFOA FREE)

NON-ABSORBABLE

MONOFILAMENT

\$99.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-PTFE3022	19mm	3/8	Silver	3-0	22"	White
S-PTFE4022	19mm	3/8	Silver	4-0	22"	White
S-PTFE5022	16mm	3/8	Silver	5-0	22"	White

POLYPROPYLENE

NON-ABSORBABLE

MONOFILAMENT

\$149.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-PP3024GB	19mm	3/8	Black	3-0	24"	Blue
S-PP4024GB	19mm	3/8	Black	4-0	24"	Blue

NYLON

NON-ABSORBABLE

MONOFILAMENT

\$149.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-NLN3024GB	19mm	3/8	Black	3-0	24"	Black
S-NLN5024GB	19mm	3/8	Black	5-0	24"	Black

SILK

NON-ABSORBABLE

BRAIDED

\$149.⁹⁵



SKU	Needle Length	Circle	Needle Color	SUT Diameter	SUT Length	SUT Color
S-SLK3024GB	19mm	3/8	Black	3-0	24"	Black
S-SLK5024GB	19mm	3/8	Black	5-0	24"	Black



It all started with a spring...

For over 170 years, dental forceps were all the same. As the dental industry continued to innovate and invent new technologies, the physical instruments clinicians use every day remained overlooked. Clinicians learned to work around hand pain, fatigue, and dangerous inadequacies. **That didn't sit right with us.**

We started with forceps.


FRINGS® introduced a patented internal spring mechanism that eliminated hand strain and improved control during extractions. It wasn't just a tweak — it was a complete overhaul of an instrument that had been overlooked for generations.

Co-founder & CEO Haseeb Sajid hit the road — going door to door, introducing FRINGS to dentists one by one. No marketing team. No distributor. Just an **unshakable belief in better dentistry.**

The hard work paid off. FRINGS caught on. And TBS Dental was born.

Designed by Dentists, For Dentists. Today at TBS Dental, we continue to disrupt alongside industry-leading clinicians. Sprouted from surgical instrumentation, our portfolio now includes streamlined solutions for critical pre- and post-extraction care.

What began as a single innovation has become a movement. One built on purpose, precision, and a relentless drive to raise the standard of care.



Our mission is simple:
Make Dentistry More Efficient.

Learn more
about TBS
Sutures

